BEFORE & AFTER RESULTS











WE BELIEVE YOUR SKIN DESERVES THE BEST

We choose industry leading German engineered microneedling technology for the best in client comfort, safety and enhanced treatment results.

TWIST treatments are ideal for:
Fine line and wrinkle reduction
Skin rejuvenation and healing
A volumised, plumper looking complexion
Acne and facial scarring reduction
Stretch mark and body scarring reduction
Increased skin elasticity and firmness
Improvements in skin texture and tone



Enhance your results with XCellaris PRO. This range of professional skin care is designed to work in complete synergy with microneedling to empower your skin at the cellular level.





Collagen Induction Therapy

Clinically proven microneedling skin rejuvenation



FREQUENTLY ASKED QUESTIONS

What happens during treatment?

Surgical-grade microneedles create tiny wounds that accelerate the body's natural healing process and work to repair the skin through increased production of collagen and elastin.

Are treatments painful?


Thanks to the optimal needle count, needle design and speed of treatment, many clients find treatments with TWIST involve minimal discomfort and are more tolerable than treatments with standard microneedling devices.

What happens after treatment?

Depending on treatment intensity you may experience redness and minor pin-point bleeding immediately after a microneedling treatment. Post-treatment redness may appear like a mild sunburn and last up to 3 days. Your clinician will recommend an aftercare plan and a post-treatment skin care product to assist in skin hydration and accelerating your results.


How many treatments are required?

Clients are able to see results in as little as 13 days with collagen production formulating 6 weeks post-treatment. Depending on the treatment area and desired outcome, most clients will require 4-6 treatments to achieve their best results.




1. During Treatment

Surgical-grade microneedles reach beyond the epidermis into the dermis, causing controlled trauma and triggering the body's natural healing response.



2. Healing Cascade

Within minutes the healing cascade starts, generating collagen-producing fibroblasts that encourage the skin to heal and renew itself from within.



3. Younger Looking Skin

New collagen and elastin forms to firm, smooth and refine the treatment area for healthier and younger looking skin. Texture and tone are improved with each treatment.

INDUSTRY LEADING TECHNOLOGY

Comfortable Treatments

Our microneedling treatments with TWIST utilise strategically designed German engineered needles that deliver up to 150 precise insertions per second which are actively retracted from the skin between insertions. This eliminates unnecessary micro-tears and damage to the skin to deliver a safer, more comfortable, and more effective treatment than when compared to other microneedling devices.

Regulated Technology

TWIST medical-grade microneedling technology is included on the Australian Register of Therapeutic Goods (ARTG #332420, #332421, #332265). Both the device and needles have met the highest European and international standards of manufacturing and certification including ISO 13485, MDSAP and CE Mark.

Consistent Results

TWIST is equipped with the most powerful microneedling handpiece on the market and unlike many other devices, is built to operate at a consistent level of power and performance that does not rely on battery power. This ensures optimal and more even treatment results in less time and with less discomfort.

Enhanced Skincare Infusions

Microneedling produces microchannels that transport our targeted corrective skincare products into the deepest layers of the skin. This allows for increased absorption and effectiveness of active skincare ingredients. The XCellaris PRO by Dermaroller skin care range is created by the makers of TWIST to work in complete synergy with microneedling treatments, empowering your skin at the cellular level.